

Application Note M15.0

Metal Scavenging using SiliaPrep[™] Thiol in Pictures

SCAVENGING of 1 000 ppm of Pd(OAc)₂ in THF USING SILIAPREP THIOL (SPE CARTRIDGE)

STEP 1: Select the appropriate SiliaPrep Thiol format and condition the cartridge.



STEP 2: Load the metal containing solution onto the SiliaPrep Thiol and let elute by gravity.



STEP 3: Metal free containing solution is coming out the SiliaPrep Thiol (wash the cartridge with additional solvent for total recovery).

